

# UCP™ Stem

## Femoral Hip System



# UCP™ Stem

## Tri-tapered Geometry

- Designed for optimal cement fixation and increased torsional stability

## Adjustment Markers

- Assist the surgeon in adjusting insertion depth of the stem to help achieve optimal leg length

## 12/14 Neck Taper

- Accommodate a complete selection of femoral heads

## Standard and High Offset

- Designed to allow optimization of soft tissue tensioning and joint stability



## Polished Surface

- Helps minimize friction at the stem-cement interface, and reduce the potential for cement failure.

## Forged Cobalt Chrome

- Designed to offer increased mechanical strength



## Centralizer

- Allows stem creeping to stable position

## Cement Restrictor



















- Prevents cement leakage

- For primary cases: Standard stem 125 mm length
- For revision cases: Long stem options in 180 mm, 210 mm, and 230 mm lengths



Size	A Medial Length	B Offset		C Vertical Height	D Neck Length		E Lateral Length
		Standard	High Offset		Standard	High Offset	
#0	110	32	37	29	41.8	48.3	132
#1	125	34	39	29	44.4	50.9	147
#2	125	36	41	29	47.0	53.5	147
#3	125	37	42	29	48.3	54.8	147
#4	125	38	43	29	49.6	56.1	147
#5	125	40	45	29	52.2	58.7	147
Long Stem							
#2	180	41		29	53.5		202
#3	180	42		29	54.8		202
#3	210	42		29	54.8		232
#3	230	42		29	54.8		252

Unit: mm

						
						
						
<b>Neck Angle</b>	130°	145°	130°	125°	131° / 125°	125°
<b>Material</b>	CoCr	CoCr	Stainless Steel	Stainless Steel	Forged CoCr / Stainless Steel	CoCr
<b>Neck Length adjustments (depth markings)</b>	15mm	15mm	10mm	8mm	11mm	8mm
<b>Tapered Design</b>	Triple	Double	Triple	Double	Triple	Triple
<b>Primary Size</b>	110, 125mm (12 sizes)	85, 95, 105, 130mm (18 sizes)	103, 106, 110.5, 115, 119.5, 124, 128.5, 133 & 137.5mm (21 sizes)	95, 115, 125 &150mm (25 sizes)	120 & 135mm (12 sizes)	127 & 149mm (9 Sizes)
<b>Revision Size</b>	180, 210 & 230mm (4 sizes)	180, 200, 230 & 260mm (7 sizes)	145.5, 148.5, 175.5 & 210.5mm (6 sizes)	205, 220, 240 & 260mm (6 sizes)	No Revision Options	No Revision Options
<b>Offset Options</b>	Standard/ High	Standard/ High/ Extra High	Standard/ High	Standard	Standard/ High	Standard

THANK YOU