



# Locking Cage™

## Revision Acetabular Hip System

# Locking Cage™

## Multiple Screw Holes

- Facilitate the surgeon's ability to utilize a desired area for screw fixation

## Pre-bent Flanges

- Designed to optimize fit with anatomy

## Modular Ischial Flange

- Provides additional support options

## Modular Hook

- Designed to provide additional prosthetic stability
- Available in 3 sizes

## Poly-axial Screw

- Designed to secure the cage to the acetabulum
- Screws are available in 15-75 mm lengths



## Locking Cage, Full XPE Cup Compatibility Guide

(Recommended maximum size, for 2.5 mm cement mantle)

Cage Size Cup Size	50 mm	54 mm	58 mm	62 mm	66 mm	70 mm
42 mm	●					
44 mm						
46 mm		●				
48 mm						
50 mm			●			
52 mm						
54 mm				●		
56 mm						
58 mm					●	
60 mm						
62 mm						●

THANK YOU